OSMOLARITY TEST CARDS

Instructions for Use

930088 REV N 2024 Trukera Medical

Trukera[™]
MEDICAL

TABLE OF CONTENTS

ENGLISH	4
БЪЛГАРСКИ	7
ČESKY	11
DANSK	15
NEDERLANDS	18
SUOMI	22
FRANÇAIS	26
DEUTSCH	30
ΕΛΛΗΝΙΚΑ	34
ITALIANO	38
日本語	42
한국어	46
NORSK	50
PORTUGUÊS	53
РУССКИЙ	57
简体中文	61
SLOVENČINA	64
ESPAÑOL	68
SVENSKA	72
TÜRKÇE	75
עברית	84
PORTUGUÊS	85
POLSKI	89

For professional in vitro diagnostic use. CLIA Complexity: Waived

CLIA CONSIDERATIONS

The TearLab® Osmolarity System is a CLIA-waived test system for human tears. Each laboratory or testing site using the TearLab Osmolarity System must have a CLIA Certificate of Waiver before starting testing. To obtain a Certificate of Waiver, call your state department of health or TearLab Customer Support at 1-858-455-6006 for an application (Form CMS-116). TearLab can provide the phone number of your state department of health and assist you in filling out the application. The TearLab Osmolarity System is a waived test so long as it is used according to the instructions in the Package Insert. Any modification by the laboratory to the test system or the test system instructions will result in this test no longer meeting the requirements for waived categorization. A modified test is considered to be high complexity and is subject to all applicable CLIA requirements.

INTENDED USE

The TearLab Osmolarity System is intended to measure the osmolarity of human tears to aid in the diagnosis of Dry Eye Disease in patients suspected of having Dry Eye Disease, in conjunction with other methods of clinical evaluation.

SUMMARY AND EXPLANATION OF THE TEST

Tears fulfill an essential role in maintaining ocular surface integrity, protecting against microbial challenge, and preserving visual acuity.¹ Hyperosmolarity has been described in the literature as a primary marker of tear film integrity.² When either the quantity or the quality of secreted tears is compromised, the total osmolarity of the tear film is increased due to increased evaporation. As a result, basal tear equilibrium is shifted to a saltier solution, which places stress on the corneal epithelia and conjunctiva. The TearLab® Osmolarity Test Card, in conjunction with the TearLab Osmolarity System, provides a quick and simple method for determining tear osmolarity using nanoliter (nL) volumes of tear fluid collected directly from the ocular surface.

PRINCIPLES OF THE PROCEDURE

The TearLab Osmolarity test utilizes a temperature-corrected impedance measurement of tear fluid to provide an indirect assessment of osmolarity.³ A voltage is applied to the tear fluid, and the electrical impedance of the dissolved tear fluid particles is monitored over time. After application of a calibration curve to the steady-state electrical impedance of the tear fluid, osmolarity is calculated and displayed as a quantitative numerical value.

MATERIALS

TearLab Osmolarity Test Card: A single-use, individually packaged, nonsterile, polycarbonate microchip containing (a) a microfluidic channel to collect 50 nanoliters (nL) of tear fluid by passive capillary action, and (b) gold electrodes embedded in the polycarbonate card to enable measurement of the impedance of the tear fluid sample in the microfluidic channel. Each Test Card is clinically hygienic and contains a protective cover. Test Cards are designed to work in conjunction with the TearLab Osmolarity System.

REAGENTS

Test Cards contain no reagents or chemicals.

Materials provided: Osmolarity Test Cards

Materials Not Provided: TearLab Osmolarity System, consisting of:

Reader, Pens, and Electronic Check Cards

Accessories: Normal and High Osmolarity Control Solutions

PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. Leave the Test Card in its sealed pouch until just before use.

3. Do not remove the protective cover until the Test Card is attached onto a Pen and immediately prior to tear collection.
4. Avoid touching the tip of the Test Card.
5. Do not use a Test Card past its expiration date.
6. Test Cards are single-use only.
7. Use only at ambient temperature of 15° - 30°C/59° - 86°F.
8. The American Academy of Ophthalmology (AAO) states that “human tears are not considered to contain significant amounts of bloodborne pathogens, and thus do not require OSHA’s [Occupational Safety and Health Administration’s] bloodborne pathogens precautions; but exposure to human tears does require good office hygiene practices such as hand washing. However, contact with tears contaminated with blood, such as in minor surgery, requires the use of bloodborne pathogen precautions.”
9. TearLab is designed to collect samples of tear fluid from the eye, a nonsterile environment. The AAO has issued guidance to minimize the transmission of ocular surface infectious agents.⁴ Prevention of transmission of these pathogens requires good hygienic techniques, such as hand washing and the cleanliness of instruments that would contact the eye. TearLab Osmolarity Test Cards are single-use, are clinically hygienic and contain a protective cover, and are individually packaged. No cleaning is required.

Although the device is manufactured to a level of disinfection deemed appropriate by the CDC (Division of Healthcare Quality Promotion; “Disinfection and Sterilization of Patient-Care Equipment, 1985”; last modified June 22, 2005), it is not sterile and there is a minimal risk of infection, ocular abrasion, or vision loss, as there is with any foreign object that is brought in contact with the eye. Appropriate hygiene and care should be adhered to.

STORAGE AND STABILITY

Store Osmolarity Test Cards at controlled room temperature (20° - 25°C/68° - 77°F); excursions are permitted to (15° - 30°C/59° - 86°F).

IMPORTANT: Test Cards not stored at room temperature should be at ambient temperature before use. Cold Test Cards may require several minutes to reach ambient temperature. Osmolarity Test Cards are stable until the expiration date marked on the label.

CALIBRATION DATA

To determine clinical performance for tear film hyperosmolarity in the diagnosis of Dry Eye Disease, a meta-analysis was performed on historical published data for tear osmolarity in samples of Normal and Dry Eye subjects.⁵ An osmolarity referent 316 mOsm/L was found to yield sensitivity of 69%, specificity of 92%, and an overall predictive accuracy of 82% for the diagnosis of Dry Eye Disease. Studies in the meta-analysis used earlier osmolarity devices, not TearLab.



PERFORMANCE OF OSMOLARITY IN META-ANALYSIS				
	Normal	Dry Eye	Total	
≤316	750	192	942	80% NPV
>316	65	429	494	87% PPV
Total	815	621	1,436	
	Specificity 92%	Sensitivity 69%		

PERFORMANCE ON PATIENTS WITH OBJECTIVE SIGNS OF DRY EYE

140 subjects were enrolled in a multicenter study (n = 45 Normal, n = 95 Dry Eye). To qualify as a Dry Eye patient, subjects were required to have a positive score on the Ocular Surface Disease Index (OSDI) and 2 or more positive indications of Tear Film Breakup Time (TBUT), Schirmer Test, Corneal Staining, Conjunctival Staining, or Meibomian Gland Dysfunction. Performance of the TearLab® Osmolarity System using these selection criteria is shown in the table below.

TEARLAB OSMOLARITY DIAGNOSTIC PERFORMANCE FOR DRY EYE DISEASE				
	Normal	Dry Eye	Total	
≤316	32	34	66	48% NPV
>316	13	61	74	82% PPV
Total	45	95	140	
	Specificity 71%	Sensitivity 64%		

Figures 1 & 2. Distribution of Osmolarities in Normal and Dry Eye Disease subjects

* Dry Eye distribution is not normal. It is skewed toward 400 mOsm/L.

OSMOLARITY WAIVER DATA

An "untrained user" study was conducted in which participants were given only the test instructions and asked to perform testing of 360 total masked, randomized samples. The samples consisted of tear matrix solution aliquoted across four levels over normal, mild, moderate and high osmolality. The participants were not given any training on the use of the test. A total of 11 participants were enrolled from three sites, representing a diverse demographic (education, age, gender, etc.) population.

RESULTS OF UNTRAINED USER STUDY			
	Trained User Average	UO Average	Regression Analysis (N=360 samples)
Level 1	295.3 ± 7.5	295.2 ± 6.6	$y = 0.98x + 6.51$
Level 2	323.3 ± 8.4	322.9 ± 7.9	$r^2 = 0.9503$
Level 3	339.6 ± 7.4	339.1 ± 6.7	Slope 95% CI = (0.95, 1.00)
Level 4	377.2 ± 7.2	377.0 ± 7.1	Intercept 95% CI = (-1.3, 14.3)

CONTRAINDICATIONS

Human tear fluid samples may be used. Collect tear fluid samples directly from the eye. Refer to the TearLab Osmolarity System User Manual for proper tear collection technique. Tear collection should be performed at the lateral (temporal) extent of the eyelid where the risk of inadvertent

injury to the cornea can be minimized, rather than adjacent to the cornea where injury is more likely.

- Do not collect tear fluid from a patient within two hours of medicinal eye drop use or use of topical medications.
- Do not collect or store tear fluid samples for transport or testing at a later time.
- Do not collect tear fluid after ocular surface staining.
- Do not collect tear fluid within 15 minutes of use of anesthetic or mydriatic (dilating) eye drops or after other invasive ocular diagnostic testing.
- Do not collect tear fluid within 15 minutes after a slit lamp examination.
- Do not collect tear fluid within 15 minutes from a patient who has been crying.

PATIENT TESTING PROCEDURE

(Refer to the TearLab Osmolarity System User Manual for more information.)

1. Turn the Reader on.
 2. Remove a Test Card from its package.
 3. Attach a Test Card by sliding the wings of the Test Card onto the TearLab Pen. The Pen will light up and beep when the card is attached properly. The green light will stay on until you collect tears or the Pen times out (after two minutes).
 4. While holding the Test Card wings, remove the protective cover from the Test Card just before tear collection.
 5. Collect a tear fluid sample as described in the TearLab Osmolarity System User Manual or Quick Reference Guide.
- IMPORTANT:** Any Test Card that does not contain a protective cover should not be used for patient testing.
6. Dock the Pen into the Reader within 40 seconds. Do NOT remove the Test Card from the Pen before docking, or all data will be lost.
 7. Locate the code on the top of the Test Card. Press the RECALL key below the UP and DOWN arrows on the Reader keypad to match the Test Card code. Press OK or wait eight seconds for the Reader to accept the code.
 8. The Reader will display the Test Result in a few seconds. Record the Test Result in the patient's chart.
 9. Remove the Test Card by pressing your thumb forward on top of the Test Card. Do not pull from the wings. Dispose in an appropriate container (refer to number eight in PRECAUTIONS).

QUALITY CONTROL

Electronic Check Card

Run the Electronic Check Card on each Pen before each day of patient testing to confirm that the system is within the manufactured calibration specifications. Refer to the Electronic Check Cards/Quality Control section in the TearLab Osmolarity System User Manual for more information.

Control Solutions

Good laboratory practice suggests the use of Normal and High Osmolality Control Solutions to ensure that 1) Test Cards are working, and 2) the test is being performed correctly.

Test both levels of control solution with each new shipment of Test Cards (even if the lot number is the same), with each new lot number, and monthly to check storage. Read the control solution instruction sheet for expected values.

- Attach a Test Card to the Pen (refer to the TearLab Quick Reference Guide).
- Do NOT collect tears.
- Instead of collecting tears, use a control solution.
- Use the blue sleeve to snap off the top of an ampule.
- Turn the ampule upside down (the fluid will not spill out).
- Touch the tip of the Pen to the control solution.
- Return the Pen to the Reader, and enter the code.
- Check control result to the expected value.
- If within the expected range, patient testing may proceed.

• If not within the expected range, you should not perform patient testing. Contact TearLab Customer Support in your area or call TearLab at (858) 455-6006.

To ensure proper performance with the TearLab Osmolarity System, only TearLab Osmolarity Control Solutions should be used. Osmolarity Control Solutions are not included with the TearLab Osmolarity System or the TearLab Osmolarity Test Cards. Osmolarity Control Solutions are single-use ampules; they should be kept at room temperature and used within their expiration date.

RESULT INTERPRETATION

TearLab test results are displayed on the LCD of the Reader in units of osmolarity or mOsm/L. No calculations are required. Osmotic concentration determinations are often expressed as either osmolarity (milliosmoles/liter, i.e., mOsm/L, or as centiosmoles/liter, i.e., cOsm/L) or osmolality (milliosmoles/kilogram, i.e., mOsm/kg, or as centiosmoles/kilogram, i.e., cOsm/kg). In tear fluid the difference between osmolarity and osmolality is less than 1%. It is common in the clinical literature to use the terms interchangeably, as the difference is insignificant.⁶ TearLab measurement range is linear from 275 mOsm/L to 400 mOsm/L. Test results outside this measurement range will be reported as either "Below Range," indicating a measurement below 275 mOsm/L, or "Above Range," indicating a measurement above 400 mOsm/L.

LIMITATIONS OF THE PROCEDURE

The results of the TearLab Osmometer should be evaluated with all clinical and laboratory data available. If the results do not agree with the clinical evaluation, additional tests should be performed. Osmolarities above or below the measurement range are very rare and should generally be confirmed with a subsequent test, as values outside the measurement range may indicate an error (e.g., user error during the test). Test only on human tears or TearLab Osmolarity Control Solutions. If no result is obtained or you encounter device problems, contact TearLab Customer Support in your area or call TearLab at (858) 455-6006.

EXPECTED RESULTS

Reference tear osmolarity values for Normal and Dry Eye Disease patients: Normal Mean 309.9 mOsm/L \pm 11.0 (288–331 mOsm/L; 90% CI 288–331)

Dry Eye Disease: Mean 324.3 mOsm/L \pm 20.1 (291–382 mOsm/L; 90% CI 284–392)

Osmolarity may differ from left to right eye, and each eye should be tested and assessed to determine which eye represents the higher osmolarity.⁷

HUMAN FACTOR AND SAFETY

Human factor usability and tear sampling discomfort and safety were evaluated at three physician office sites. Results concluded that the TearLab Osmolarity System has a positive usability factor of 4.8 out of 5.0. Among 234 subjects, the average discomfort index (1=none to 4=severe) was 1.14. Physician examination pre- and post-tear collection indicated no trauma to the cornea or conjunctiva.


PRECISION The results below come from a TearLab instrument calibrated against a reference standard solution prepared from dried, high-purity sodium chloride traceable to the National Institute of Standards and Technology (NIST).

AVERAGE OSMOLARITY (MOSMS/L)	WITHIN RUN (SD)	WITHIN RUN (CV%)	TOTAL (SD)	TOTAL (CV%)
280	3.8	1.34%	5.2	1.87%
294	5.5	1.85%	7.3	2.47%
316	4.5	1.41%	6.6	2.08%
345	4.5	1.30%	8.0	2.33%

ACCURACY (METHOD COMPARISON)


The correlation study was performed internally by the manufacturer using contrived tear samples of various osmolarity levels within the clinical reference range on both the TearLab Osmolarity System and the Wescor Model 5520 vapor pressure osmometer calibrated to NIST-traceable standards.

NO. SITES	N	REGRESSION LINE	R ²
1	80	$y = 0.9146x + 23.061$	0.9443

At each of the three physician office sites, 40 contrived tear specimens across seven levels of the clinically significant range were prepared and measured on the TearLab Osmolarity System. The physician office laboratories did not have access to the Wescor 5520 Vapro[®] vapor pressure osmometer. Wescor values were determined from an average of two to three measurements on each level of osmolarity immediately prior to the beginning of the study.

NO. SITES	N	REGRESSION LINE	R ²
3	120	$y = 0.9402x + 12.512$	0.9515

INTERFERING SUBSTANCES

Tear samples are to be collected directly from the eye. Do not collect tears within two hours of medicinal or staining eye drop use or within 15 minutes of use of anesthetic or mydriatic (dilating) eye drops. Do not collect tears within 15 minutes from patients who have been crying.

FDA MEDWATCH

Report a serious adverse event, product quality problem, product use error, or therapeutic inequivalence/failure that you suspect is associated with the use of the TearLab Osmolarity System to TearLab Customer Support (Tel: 858-455-6006) and/or to FDA MedWatch (Tel: 800-FDA-1088), (Fax: 800-FDA-0178), or (www.fda.gov/medwatch).

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This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

This page was intentionally left blank for your notes.

For ordering information, contact your local sales representative or TearLab Customer Support in your area.

1. Sullivan DA, Dartt DA, et al. Lacrimal Gland, Tear Film, and Dry Eye Syndromes 2: Basic Science and Clinical Relevance. *Advances in Experimental Medicine and Biology*. 1998; 438.
2. Definition and Classification of Dry Eye. Report of the Diagnosis and Classification Subcommittee of the Dry Eye Workshop (DEWS). *The Ocular Surface*. 2007;5(2):75-92.
3. Pedersen-Bjergaard K, Smidt BC. Electrolytic Conductivity, Osmotic Pressure, and Hydrogen Ion Concentration of Human Lacrymal Fluid. *Acta Derm Venereol Suppl* (Stockh). 1952;32(29):261-7.
4. Minimizing Transmission of Bloodborne Pathogens and Surface Infectious Agents in Ophthalmic Offices and Operating Rooms. American Academy of Ophthalmology, Information Statement, Available at: <http://one.aao.org/CE/PracticeGuidelines/ClinicalStatements.aspx?p=3> Accessed March 1, 2003.
5. Tomlinson A, Khanal S, Ramaesh K, Diaper C, McFadyen A. Tear Film Osmolarity: Determination of a Referent for Dry Eye Diagnosis. *Investigative Ophthalmology & Visual Science*. October 2006;47(10): 4309-4315.
6. Erstad BL. Osmolality and Osmolarity: Narrowing the Terminology Gap. *Pharmacotherapy*. 2003;23(9):1085-6.
7. Lemp MA, Bron AJ, Baudouin C, et al. Tear Osmolarity in the Diagnosis and Management of Dry Eye Disease. *Am J Ophthalmol* 2011: ePub February 2011.

 Trukera Medical, 142309 Winchester Rd. Ste. I, Temecula, CA 92590 USA
1-858-455-6006, www.trukera.com

 Emergo Europe, Westervoortsedijk 60 6827 AT Arnhem The Netherlands

 For complete instructions for use, please refer to the TearLab or ScoutPro Osmolarity System User Manual.

 For In Vitro Diagnostic Use

CE

US Pat. # 7,017,394, 7,051,569, 7,111,502, 7,129,717, 7,204,122, 8,020,433, 7,987,702, 7,905,134, 7,810,380, 7,574,902

930088 REV N 2024 Trukera Medical

Trukera[™]
M E D I C A L